

Standard	Dilution	Contact time	Temperature	Strains
EN 1276 - Bactericidal	ready to use	5 minutes	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae, Escherichia coli</i>
EN 1499 - Bactericidal	ready to use	1 minute	20°C	<i>Escherichia coli</i>
EN 13727 - Bactericidal	ready to use	1 minute	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus, Escherichia coli K12, Enterococcus hirae</i>
EN 13624 - Yeasticidal	ready to use	1 minute	20°C	<i>Candida albicans</i>

Material to be used exclusively by the Sutter Professional sales force