

Standard	Dilution	Contact time	Temperature	Strains
EN 1040 - Bactericidal	3%	5 minutes	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus</i>
EN 1276 - Bactericidal	3%	5 minutes	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae, Escherichia Coli</i>
EN 1276 - Additional Bactericidal	3%	5 minutes	20°C	<i>Salmonella enterica subs, Listeria monocytogenes</i>
EN 13697 - Additional Bactericidal	3%	15 minutes	40°C	<i>Salmonella typhimurium, Listeria monocytogenes</i>

Material to be used exclusively by the Sutter Professional sales force