

Standard	Dilution	Contact time	Temperature	Strains
EN 1276 - Bactericidal	ready to use	5 minutes	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae, Escherichia Coli</i>
EN 1650 - Yeasticidal	ready to use	5 minutes	20°C	<i>Candida Albicans</i>
EN 13697 - Bactericidal	ready to use	5 minutes	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae, Escherichia Coli</i>
EN 13697 - Yeasticidal	ready to use	5 minutes	20°C	<i>Candida Albicans</i>
EN 14476+A1 - Virucidal	ready to use	5 minutes	20°C	<i>Virus A/H1N1</i>
EN 1276 - Additional Bactericidal	ready to use	5 minutes	20°C	<i>Salmonella enterica subsp. salamae serotyoe Tranoroa, Listeria monocytogenes</i>
EN 13697 - Additional Bactericidal	ready to use	5 minutes	40°C	<i>Salmonella typhimurium, Listeria monocytogenes</i>
EN 13727 - Bactericidal	ready to use	5 minutes	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae</i>
EN 13624 - Yeasticidal	ready to use	15 minutes	20°C	<i>Candida Albicans</i>

Material to be used exclusively by the Sutter Professional sales force