

Standard	Dilution	Contact time	Temperature	Strains
EN 1276 - Bactericidal	ready to use	5 minutes	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae, Escherichia Coli</i>
EN 1650 - Yeasticidal	ready to use	15 minutes	20°C	<i>Candida Albicans</i>
EN 1650 - Additional Fungicidal	ready to use	15 minutes	20°C	<i>Penicillum Aurantiogriseum</i>
EN 13697 - Bactericidal	ready to use	5 minutes	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae, Escherichia Coli</i>
EN 13697 - Yeasticidal	ready to use	15 minutes	20°C	<i>Candida Albicans</i>
EN 13697 - Additional Fungicidal	ready to use	15 minutes	20°C	<i>Penicillum Aurantiogriseum</i>
EN 13727 - Bactericidal	ready to use	5 minutes	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae</i>
EN 13624 - Yeasticidal	ready to use	15 minutes	20°C	<i>Candida Albicans</i>

Material to be used exclusively by the Sutter Professional sales force