

Standard	Dilution	Contact time	Temperature	Strains
EN 1040 - Bactericidal	2-4%	5 minutes	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus</i>
EN 1276 - Bactericidal	2-4%	5 minutes	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae, Escherichia Coli</i>
EN 13697 - Bactericidal	2-4%	5 minutes	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae, Escherichia Coli</i>
EN 13727 - Bactericidal	5%	5 minutes	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus, Escherichia coli K12, Enterococcus hirae</i>
EN 13624 - Yeasticidal	5%	15 minutes	20°C	<i>Candida albicans</i>

Material to be used exclusively by the Sutter Professional sales force