

Standard	Dilution	Contact time	Temperature	Strains
EN 1276 - Bactericidal	3%	5 minutes	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae, Escherichia coli</i>
EN 1650 - Fungicidal	3%	15 minutes	20°C	<i>Candida albicans, Aspergillus niger</i>
EN 13697 - Bactericidal	3%	5 minutes	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae, Escherichia coli</i>
EN 13697 - Fungicidal	3%	15 minutes	20°C	<i>Candida albicans, Aspergillus niger</i>
EN 1276 - Additional Bactericidal	5%	5 minutes	20°C	<i>Salmonella enterica subsp. salamae serotyoe Tranoroa, Listeria monocytogenes</i>
EN 1276 - Additional Bactericidal	4%	20 minutes	20°C	<i>Clostridium Difficile</i>
EN 13697 - Additional Bactericidal	5%	5 minutes	40°C	<i>Salmonella typhimurium, Listeria monocytogenes</i>
EN 13697 - Additional Bactericidal	4%	20 minutes	20°C	<i>Clostridium Difficile</i>
EN 14476 - Virucidal	3%	5 minutes	20°C	<i>Murine Norovirus, Poliovirus Type 1, Adenovirus Type 5</i>
EN 13624 - Yeasticidal	3%	15 minutes	40°C	<i>Candida albicans</i>
EN 13727 - Bactericidal	10%	5 minutes	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus, Escherichia coli K12, Enterococcus hirae</i>

Material to be used exclusively by the Sutter Professional sales force