

Standard	Dilution	Contact time	Temperature	Strains
EN 13697 - Bactericidal	6%	5 minutes	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae, Escherichia Coli</i>
EN 13697 - Fungicidal	1,5%	15 minutes	20°C	<i>Candida Albican, Aspergillus brasiliensis</i>
EN 13727 - Bactericidal	6%	5 minutes	30°C	<i>Staphylococcus Aureus, Pseudomonas Aeruginosa, Enterococcus Hirae</i>
EN 13624 - Yeasticidal	1,5%	15 minutes	30°C	<i>Candida Albicans</i>
EN 14476 - Virucidal	6%	5 minutes	20°C	<i>Murine Norovirus, Poliovirus Type 1, Adenovirus Type 5</i>

Material to be used exclusively by the Sutter Professional sales force