



TEST AND CONDITIONS

XTRA-CLOR

Standard	Dilution	Contact time	Temperature	Strains
EN 13697 - Bactericidal	ready to use	5 minutes	20°C	<i>Pseudomonas Aeruginosa, Staphylococcus Aureus, Enterococcus Hirae, Escherichia Coli</i>
EN 13697 - Fungicidal	ready to use	15 minutes	20°C	<i>Candida Albicans, Aspergillus Brasiliensis, Penicillium Chrysogenum</i>

Material to be used exclusively by the Sutter Professional sales force

